

	<b>HOLISTA TRANZWORLD PRIVATE LIMITED</b> 2/91, MARAVANMADAM, ANTONIYARPURAM, TUTICORIN -628 101	Doc No : HTPL-SOP23
	<b>Food Safety and Quality Management System</b>	Issue/Rev : 1.0
	<b>Lab Test Procedure SOP</b>	Date : 19.09.2022

## 9.DETECTION OF *SALMONELLA*

### SCOPE:

To detect *Salmonella* in Coconut milk and Desiccated Coconut products.

### PRINCIPLE:

*Salmonella* are facultative anaerobic, Gram-negative, non-spore-forming, rod shaped and predominantly motile enterobacteria. Measured aliquots of samples are pre-enriched in non-selective enrichment medium (BPW) and incubated at  $37 \pm 1^\circ\text{C}$  for 16 h to 20 h. After incubation measured aliquots are inoculated into selective enrichment medium (RV medium & SC medium) and incubated RV medium at  $42 \pm 1^\circ\text{C}$  for 24 h and SC broth at  $37 \pm 1^\circ\text{C}$  for 24 h. After incubation sub cultures are made on BGA and XLD and incubated at  $37 \pm 1^\circ\text{C}$  for 24 h to 48 h. The presence of *Salmonella* in the sample is determined by observing the characteristic colonies developed on the plates.

### DILUENT, CULTURE MEDIA AND REAGENTS :

- Buffered peptone water (BPW)
- Rappaport-Vassiliadis medium (RV)
- Selenite Cystine medium (SC)
- Brilliant green agar (BGA)
- Xylose lysine deoxycholate agar (XLD)

### APPARATUS, INSTRUMENTS AND GLASSWARES:

- Laminar airflow chamber
- Autoclave
- Hot air oven
- Sterile petri plates
- Weighing balance
- Dilution bottles
- Water bath
- Screw cap tubes
- Sterile spatula
- Micropipette
- Sterile tips

	<b>HOLISTA TRANZWORLD PRIVATE LIMITED</b> 2/91, MARAVANMADAM, ANTONIYARPURAM, TUTICORIN -628 101	Doc No : HTPL-SOP23
	<b>Food Safety and Quality Management System</b>	Issue/Rev : 1.0
	<b>Lab Test Procedure SOP</b>	Date : 19.09.2022

## PREPARATION OF TEST SAMPLE:

### 1.Desiccated coconut

All the samples are aseptically mixed in order to make the sample uniform.

### 2. Coconut milk

All the samples are aseptically mixed in order to make the sample uniform

## PROCEDURE:

- Aseptically transfer 25g of sample into 225 ml of sterile BPW in screw cap bottles (For specific customer requirements increase the sample quantity with BPW in the ratio of 1: 9 e.g. 375g sample into 3375ml sterile BPW). Mix the sample by shaking 25 times in 7 seconds with a 1 foot(30cm) movement. For swab samples add 10 ml swab samples to 90 ml of sterile BPW and mix well.
- Incubate the bottles at  $37 \pm 1^{\circ}\text{C}$  for 16 to 20 h.
- After incubation transfer 10ml sample to 100ml SC medium and 0.1ml sample to 10ml RV medium.
- Incubate SC medium at  $37 \pm 1^{\circ}\text{C}$  for 24 h and RV medium at  $42 \pm 1^{\circ}\text{C}$  for 24 h.
- After incubation streak out a loopful sample from SC medium and RV medium on to BGA and XLD agar.
- Incubate the plates at  $37 \pm 1^{\circ}\text{C}$  for 24 to 48 h.
- After incubation observe the plates for the presence of *Salmonella* colonies (Pinkish white colonies on BGA and pink colonies with black center on XLD).
- If there is any typical colonies, sub culture on Nutrient agar plates and do biochemical confirmation using HI Salmonella identification kit.

## SEROLOGICAL CONFIRMATION:

- Place three separate loopfuls of normal saline (0.85% sodium chloride) on a clean glass slide.
- Take a small part of a suspect *Salmonella* colony from an overnight nutrient agar plate and mix thoroughly with both drops of normal saline on the slide to obtain a smooth suspension.
- Add one drop of Salmonella Poly O & Vi antisera to one bacterial suspension

	<b>HOLISTA TRANZWORLD PRIVATE LIMITED</b> 2/91, MARAVANMADAM, ANTONIYARPURAM, TUTICORIN -628 101	Doc No : HTPL-SOP23
	<b>Food Safety and Quality Management System</b>	Issue/Rev : 1.0
	<b>Lab Test Procedure SOP</b>	Date : 19.09.2022

and

Salmonella Poly H antisera to second bacterial suspension drops on the slide, to the other (control) add one loopful of normal saline.

- Mix the antiserum with the bacterial suspension using a sterile loop.
- Gently tilt the slide back and forth for one minute and observe for agglutination under normal lighting conditions.

#### **EXPRESSION OF RESULT:**

Report result as Present/Absent of *Salmonella*/tested volume of the sample.

#### **REFERENCE:**

IS 5887 (Part 3) : 1999 Reaffirmed 2009 / ISO 6579 : 1993