

	HOLISTA TRANZWORLD PRIVATE LIMITED 2/91, MARAVANMADAM, ANTONIYARPURAM, TUTICORIN -628 101	Doc No : HTPL-SOP15
	Food Safety and Quality Management System	Issue/Rev : 1.0
	Lab Test Procedure SOP	Date : 19.09.2022

15. Caustic Strength determination

1. Purpose:

To determine the concentration of sodium hydroxide (caustic soda) in cleaning solutions used in food processing equipment.

2. Scope:

Applicable to all cleaning solutions prepared and used in the plant using sodium hydroxide (NaOH).

4. Principle:

Sodium hydroxide reacts with a known concentration of standard hydrochloric acid in the presence of phenolphthalein indicator. The amount of HCl used allows calculation of NaOH concentration.

5. Reagents and Chemicals:

- Standard Hydrochloric acid (HCl), approx. 0.1 N
- Phenolphthalein indicator (1% in ethanol)
- Distilled water

6. Apparatus:

- Burette
- Pipette
- Conical flask
- Beaker
- Glass rod

7. Procedure:

1. Sample Preparation:

Take 10 ml of the caustic cleaning solution using a pipette and transfer into a clean conical flask. Dilute with 50 ml of distilled water.

2. Add Indicator:

Add 2-3 drops of phenolphthalein indicator. The solution will turn pink if caustic is present.

3. Titration:

Fill burette with 0.1 N HCl. Titrate the caustic solution until the pink color just disappears. Note the volume (V) of HCl used.

	HOLISTA TRANZWORLD PRIVATE LIMITED 2/91, MARAVANMADAM, ANTONIYARPURAM, TUTICORIN -628 101	Doc No : HTPL-SOP15
	Food Safety and Quality Management System	Issue/Rev : 1.0
	Lab Test Procedure SOP	Date : 19.09.2022

8. Calculation:

Normality of NaOH (N) = $(V \times N_{\text{HCl}}) / 10$

Strength (% w/v) = Normality \times 4.0

Where:

- V = Volume (ml) of HCl used
- N_{HCl} = Normality of HCl (usually 0.1 N)
- 10 = Volume (ml) of sample taken
- 4.0 = Equivalent weight factor to convert N to % for NaOH